

CELLiST_{TM} Supplement Cys2

Model	Package type	Qty.
CELLiST _{TM} Supplement Cys2	Aluminum Pouch Supplement for 1L feed medium preparation	17.0 g

Properties

- Chemically-Defined.
- Does not contain animal-derived components.
- Contains no hydrolysates, extracts or other undefined components.
- Does not contain protein like growth factors.
- Does not contain L-glutamine sources, sodium bicarbonate, or poloxamer.

<u>Note</u>: Cys2 is meant to be used as a substitute for the second component (high pH component) in 2-component feed products (for example, replacing CB7b component).

Instructions for 1L feed media preparation

- 1. Prepare proper container (beaker or flask) and stirrer. To ensure efficient stirring, a 2 L or 3 L container is recommended.
- 2. Add cell culture grade water (purified water) to the container up to 800 ml.
- 3. Add the main feed medium component in the required amount needed for preparation of 1000 mL feed.
- 4. Add <u>17.0 g of CELLiST_{TM} Supplement Cys2</u>.
- 5. Mix for at least 30 minutes. (Note: It will not completely dissolve at this stage).
- 6. Slowly add 5N NaOH in order to adjust the pH to the range of 6.6~6.8.
- 7. Mix for at least 60 minutes until completely dissolve.
- 8. Adjust the pH using either a 5N or 10N NaOH solution or HCl solution to the range of 6.6~6.8.
- 9. Mix for another 60 minutes.
- 10. Add cell culture grade water to make up final volume (1000mL) and stir for about 10 minutes.
- 11. Make sure the pH is in the correct range (6.6~6.8).
- Sterile filter in a clean bench, using a membrane filter with pore size of 0.2 to 0.22 μm in diameter (using a positive pressure system).
- 13. Keep the prepared medium refrigerated (2~8°C) until use. It is recommended to store the feed medium in 50 mL sealed Falcon tubes leaving a minimum amount of free air space (in order to minimize oxidation of the medium).

Storage

Store under cool (2~8°C), dark and dry conditions until use.

Use

For laboratory use only. Not intended for human or therapeutic use.