

Cell culture media For better lives

CELLISTTM BASAL10 / BASAL10P

Shortcut to Official Channel

CELLISTTM BASAL10 / BASAL10P

Overview

CELLiST™ BASAL10 is a high-performance growth medium specifically designed to support robust cell growth and viability in various CHO cell lines. This chemically defined medium, free from animal-derived components, ensures consistent and reliable results across different production batches. Leveraging Ajinomoto's extensive expertise in developing and manufacturing amino acids, CELLiST™ BASAL10 provides everything your CHO cell line needs for stable, high-yield protein production. Suitable for use with any CHO cell line, CELLiST™ BASAL10 offers a comprehensive solution for biotherapeutic production with superior quality and reproducibility.



| CELLIST™ BASAL10 Features | Benefits |
|--|--|
| Enriched with amino acids for enhanced cell growth. | Rapid cell proliferation, achieving higher cell densities faster. |
| Versatility across cell lines. | Compatible with a wide range CHO cell line, including CHO-K1, CHO-GS, CHO-S and their derivatives. |
| Chemically defined, protein-free medium with no animal-derived components. | Reduces risk of viral contamination, ensures batch consistency. |
| Supplied in fine powder form. | Easy to dissolve and allows for prolonged shelf life and ease of transportation. |
| Available in test samples and bulk sizes. | Flexibility according to usage requirements. |
| Suitable for batch, fed-batch, and perfusion cell cultures, at any scale. | Versatility across cell culture processes and scales. |

BASAL10 Product List

| Product | Catalog No. | Specification | Amount | Packaging Information |
|----------|-------------|--------------------|----------|-----------------------|
| BASAL10 | BS1 | - | 24.3 g/L | 1L / 50L |
| BASAL10P | BS4 | Includes Poloxamer | 25.3 g/L | 1L / 50L |

| Cell Culture Performance

Comparative studies show superior culture growth and productivity with the CELLiST™ BASAL10 + CELLiST™ F7 media combination. In a fed-batch process using a CHO-K1 cell line expressing IgG1 antibody, our media outperforms top global competitors in viable cell density and IgG titer.

Cell Growth

-CELLIST_{TM} BASAL10 + CELLIST_{TM} F7 -O-Supplier A -O-Supplier A2 -O-Supplier B -O-Supplier C

Protein Production




Figure 1. Fed-batch process was performed in a 125 mL Erlenmeyer flask system, using CHO-K1 cell line expressing IgG1 antibody. Feeding in CELLiST™ group was performed at a concentration of 6% (v/v) on day 4, 6, 8, 10 and 12.

A, A2, B and C represent basal/feed media combinations from major media suppliers. Culture and feeding manners were performed according to each manufacturer's recommendations.

Other Products

| Products | Туре | Features |
|-------------------------|------------------------|---|
| CELLIST™ F7 | Feed medium | High-performance, all-in-one feed media that easily dissolves and remains stable at neutral pH, thanks to Ajinomoto Group's unique cysteine-stabilization technology. This enhances cysteine availability, boosting performance and enabling hassle-free, single-agent feeding. |
| AminoSupplement Cys1 | Feed medium supplement | Performance-enhancing supplement for any single-agent feed media, improving culture performance and stability with Ajinomoto Group's unique cysteine stabilization technology. |
| AminoSupplement Cys2 | Feed medium supplement | Performance-enhancing supplement that eliminates the need for high-pH feeding in <u>dual-agent feed media</u> . Easily added at neutral pH, ensuring high stability and simplified processes. |
| Supplement Gly-L-Tyr | Feed medium supplement | Highly soluble tyrosine supplement that boosts tyrosine concentrations for any neutral pH feeding process, enhancing culture growth without solubility concerns. |

| Sales & Technical Support

CELLiST™ Sales Team

Tel: +82-032-210-2600 Mail: sales@ajinomotocellist.com

CELLIST™ Solution Center India (CSC-IN)

2nd Floor, Block 05, TIMES CORPORATE PARK Thaltej, Ahmedabad, Gujarat 380058, India

Tel: +91-87359-10813 Mail: parshadg.hirapara.ut3@asv.ajnomoto.com Tel: +86-21-6219-2705 Mail: mingqi_fan@ajnomoto.com.cn

CELLIST™ Solution Center Korea

70, Songdogwahak-ro, Yeonsu-gu, Incheon, Republic of Korea Mail: yaron.silberberg.dk2@asv.ajinomoto.com

CELLIST™ Solution Center China (CSC-CN)

Room 501, Building 9, NO.185 Jiangtian East Road, Songjiang District, Shanghai, China



Documents

Eat Well, Live Well.

